



TWELFTH EDITION *Catalog of Services - Update # 1*

November 2009

Test Code	Test Name	Page(s)	Update	New Information
3402	Chromosome Analysis and Fragile X Syndrome	44 and 67	CPTs, Fee	<p>CPTs: 88230, 88262, 83891, 83892 x 3, 83893, 83894, 83896, 83897, 83900, 83909, 83912</p> <p>Fee: \$1884.72 (effective November 1, 2009) <i>(pgs 6 & 12 of Fee Supplement)</i></p>
8147	Cytomegalovirus (CMV) by Rapid PCR	57	Specimen Requirement	<p>Submit only ONE of the following specimens: 5.0 mL EDTA or ACD-B whole blood (1.5 mL minimum); 1.0 mL EDTA or ACD-B plasma (0.5 mL minimum); 1.0 mL serum (0.5 mL minimum); 1.0 mL CSF, amniotic fluid or urine (0.5 mL minimum); 3 mm³ fresh tissue frozen; 1.0 mL bronchial brush/wash. Blood collected with heparin anticoagulant is NOT acceptable. Submit specimen in sterile leak proof container. INDICATE SOURCE.</p> <p>Stability: RMT – 48 hours REFT – 8 days</p> <p>NOTE: Do NOT spin CSF specimens.</p>
7542 Test Code discontinued; Use new code 7146	1,25-Dihydroxyvitamin D	58 and 112	<p>Code 7542 no longer orderable; use Code 7146 for this test</p> <p>Specimen requirement</p> <p>Methodology</p> <p>Stability</p>	<p>2.0 mL serum (1.1 mL minimum), room temperature.</p> <p>Methodology: Liquid Chromatography/Tandem Mass Spectrometry <i>Note: Reports Total, D2 and D3 of 1,25 Dihydroxyvitamin D</i></p> <p>Stability: RMT – 14 days REFT – 14 days FROZEN – 14 days</p>

2356	Echovirus (4,7,9,11,30) Antibodies	63	Specimen Requirement and Stability	1.0 mL serum only (0.5 mL minimum). Send room temperature or refrigerated. Stability: RMT – 7 days REFT – 14 days
7813	Insulin Antibodies	80	Specimen requirements Stability	Do not use serum separator or gel barrier tubes. Send serum refrigerated. Stability: RMT – 3 days REFT – 14 days
4980	Hepatitis C Virus RNA, Quantitative	74	Specimen Requirement	3.0 mL FROZEN EDTA plasma (1.1 mL minimum). Samples should be spun within 6 hours of collection. Separate plasma and freeze. Serum or ACD plasma samples are also acceptable. Sodium heparin plasma is NOT acceptable. Patient should NOT be on heparin therapy.
1950	Propafenone (Rythmol)	94 and 101	Specimen requirements	Plasma is no longer acceptable.